MATERIALI

Resine Termoelastiche

ASAKO (polipropilene)

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GB

FIELD OF APPLICATION

The thermoplastic resin **Asako** is intended for use by skilled personnel at the prosthodontist lab for the creation of partial dentures and hooks. **Asako** is a Class I product (maximum duration in the oral cavity 30 days), pursuant to Directive 93/42/EEC as amended, by applying Rule 5 in Annex IX.

MAIN FEATURES

- Easy processing
- **Asako** is compatible with all other injection systems
- Excellent fluidity during injection
- Natural appearance
- Lightness and flexibility
- Chromatic unalterability
- Stain-resistant

INSTRUCTIONS FOR USE

The material is in granules, so you can use just the amount needed.

1. Fill the cartridge with the desired quantity. It is preferable not to completely fill the cartridge so as to allow the escape of any air (leave about 1 cm of space between product and cap). Avoid direct contact with hands to avoid contamination.
2. Close the cartridge with the cap.
3. Lubricate the cartridge with silicone spray suitable for high temperatures (Universal Spray).
4. Insert the cartridge into the heat press according to the manufacturer's instructions. Set the melting temperature to 230°C and keep it at that setpoint for 12'30". Leave the pressure of use set by the heat press manufacturer, then start the cycle.
5. At the end of the cycle, wait for the muffle to completely cool without creating thermal shocks
6. Remove from the muffle, finish and polish, making sure the material does not overheat.
7. During the finishing phase, we recommend using Kit Set Z and Smile Set combined with Lux polishing.
8. Avoid treating the prosthetic product with heat sources like steam or flames.
9. For optimal injection, we recommend the wax pattern or plate to have a thickness between (no more than) 2 mm and (no less than) 1.5 mm.
10. Use 5 mm thick wax profiles for casting posts
11. We recommend using Techifit Light gypsum-resin insulation.

- Amount of material recommended according to the prosthetic product:
- 6/8 g hook or single-element spacer
 - 12 g multi-element spacer or single bridge, up to 4 elements
 - 16/18 g bilateral sample, as Class 1 (Kennedy classification)
 - 22/24 g near-total or total prosthesis

STORAGE CONDITIONS

Store the product, hermetically sealed, in a cool, dry place. Expiration guidelines for our products refer to sealed packaging. Once opened, the product is guaranteed for one year.

DISPOSAL

Dispose of according to current legislation.

WARNINGS AND CONTRAINDICATIONS

The thermoplastic resin **Asako** is a non-cytotoxic material and does not call for any special precautions during use and storage. To clean partial dentures, we recommend cleaning every day using a soft bristle brush. Do not use abrasives such as toothpaste, etc. The product is not resistant to chemical agents with oxidizing effects and strong acids. During the finishing and polishing phase, be careful not to overheat the piece. Since it is made of plastic, it could undergo changes and irritating gases could develop. The material is not visible to x-rays.

LIABILITY ON THE PRODUCT: Our consulting and the instructions are based on current knowledge and experience and should only be considered as an indication without any binding value. The user shall not be exempt from verifying that the product is suitable for the expected application and result. Any liability for damages will be limited to the value of the merchandise supplied and used by the user. Our products are constantly being updated. We reserve the right to make changes. We recommend reading the instructions carefully before using a new package.

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