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CHEMICAL-PHYSICAL AND MECHANICAL CHARACTERISTICS OF TITANIUM

Titanium is a metallic non magnetic element which is present in nature as Oxide (TiO₂).

Pure titanium (sponge) is obtained:

- 1) After TiO₂ chlorination by thermochemical reduction of TiCl₄ (Kroll and Hunter process).
- 2) After TiO₂ chlorination by electrolytic reduction of TiCl₄ into mineral salts.

The obtained titanium sponge can be remelted in blooms, ingots and billets.

By rolling, drawing and extrusion you obtain the main semimanufactures available on the market : bars, tapes, wires, sheets and pipes.

Main chemical-physical characteristics of titanium and its alloys:

- 1) High corrosion resistance
- 2) Good mechanical resistance (specially in alloys)
- 3) Lightness (4,5 kg/dm³)
- 4) Low thermal expansion coefficient
- 5) Low electric conductivity
- 6) High elasticity
- 7) Non magnetic element

Those extraordinary characteristics have increased the use of titanium for many applications.

Titanium is available on the market as pure titanium or as alloy. Pure titanium has a crystalline structure of alpha type (compact hexagonal): it is graded from gr. 1 to gr. 4 and presents an increasing mechanical resistance.

Titanium alloys have a crystalline structure of alpha+beta type (body-centered cubic), with addition of stabilizing agents, such as Aluminium and Vanadium. Those alloys are the most commonly used, since they give best ratio between weight and mechanical resistance.

The actually most commonly used alloy is Ti6Al4V (gr.5) , largely employed in screws and mechanical components manufacturing.

The following tables show:

- Chemical composition of various types of pure titanium and some of its most commonly used alloys
- Relative physical and mechanical characteristic.

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PHYSICAL AND MECHANICAL CHARACTERISTICS OF TITANIUM AND ALLOYS

GRADO/GRADE	Titanio commerciale puro Commercially pure Titanium				Leghe di Titanio alpha+beta Half+Beta Titanium Alloy		
	Grado 1	Grado 2	Grado 4	Grado 7	Grado 5	Grado 5 ELI	Grado 9
Rm (MPa)	240	345	550	375	900	900	620
Rp 0.2 (MPa)	170	275	485	345	830	830	520
elongation %	24	20	15	20	10	10	15
KCU resilience (J)	110-160	40-80	14-27	40-80	15-20	24	54
hardness	120 HB/30	160 HB/30	250 HB/30	160 HB/30	36 HRC	35 HRC	20 HRC
own weight Kg/dm3	4.51	4.51	4.51	4.51	4.40	4.40	4.48
modulus of elasticity (GPa)	103	103	107	107	110	110	103
shear modulus (GPa)	40	40	40	40	40	40	40
melting point (°C)	1668	1668			1650		
coefficient of thermal exposure (1/°C)	8.4x10 ⁻⁶	8.4x10 ⁻⁶			8.5x10 ⁻⁶		
specific heat (cal/g-°C)	0.124	0.124			0.138		
thermal conductivity (cal/cm ² -sec-°C/cm)	0.052	0.052			0.024		
electrical resistance (μ Ω-cm)	55	55			175		
weldability	excellent	excellent	excellent	excellent	good	good	excellent
forging temper. (°C)	870-930	870-930	900-930	900-930	950-980	950-980	870-900
molding temper. (°C)	815-870	815-870	840-870	815-870	900-950	900-950	760-815
Beta Transus (°C+/-25°)	882	900	940	900	1000	990	795
type of semimanufacture	Barre billetta lamiera filo tubo vergella bars, billets, sheets, wire, tubes, wire rods						
applications	In chemical and nautical industry for corrosion resistance. In aeronautical components where the highest mouldability is required.		Aircraft engine components. Fasteners, aeronautical structures, it's the most common alloy.		High tenacity to rupture, low temperatures. Commonly used for medical and marine applications.		hydraulic pipes in aeronautics. Tensile strength can be increased to 700 Mpa with cold forming

PRINCIPALI PROPRIETA' MECCANICHE DELLE VITI
MAIN MECHANICAL PROPERTIES OF SCREWS

MATERIALI	ANALISI CHIMICA	DENSITA'	CARICO ROTTURA	LIMITE SNERV.	ALLUNGAMENTO
MATERIALS	CHEM. ANALYSIS	DENSITY	TENSILE STRENGTH	YIELD POINT	LENGTHENING
		KG/dm ³	Min/max Mpa (1)	0,2% Mpa (1)	% min
TITANIUM	ASTM 348 Gr.2	4,5	345/450	275	20
TITANIUM	ASTM 348 Gr.5	4,5	895/1100	825	10
ALUMINIUM	ASTM 7075	2,7	550/650	490	12
MARAGING	Ni-Mo-Cr-Co	8	1850/2000	1750/1900	08/10/11
ALLOY	Ni-Co-Cr-Mo	10	2000	1950	8

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POSSIBILITY OF SUBSTITUTION OF THE STEEL SCREWS WITH OUR PRODUCTS

		3.6	4.6	4.8	5.6	5.8	6.8	8.8		10.9	12.9
								<=M16	>=M16		
Unitary tensile Rm N/mm ²	nom.	300	400	400	500	500	600	800	800	1000	1200
	min.	330	400	420	500	520	600	800	830	1040	1220
Unitary yeld point Rel N/mm ²	nom.	180	240	320	300	400	480				
	min.	190	240	340	300	420	480				
Unit. Dev. Of prop. Rp02 N/mm ²	nom.							640	640	900	1080
	min.							640	640	940	1100
Vickers hardness	nom.	95	120	130	155	160	190	250	258	320	375
HV.F>=98 N	max.	220	220	220	220	220	250	310	320	375	430
Sostituzione con Tit. GR5		X	X	X	X	X	X	X	X	X	
Sostituzione con AI 7075		X	X	X	X						

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Biotan™ Nb Alloy

Titanium Modellation

The modellation is done in the conventional method, as described in our **titanium information brochure** and the **sprueing tips**.

Investing

Biotan™ Nb must only be invested with the **special investment material Biotan MG** (for frameworks) or **Biotan Vest MG** (for frameworks). The according heating rates must be kept as listed in the instructions for these investment materials.

Preheating Biotan MG/Vest MG

Temperature increase 5°C per minute after each holding point

1. holding point at **300°C** holding time: **45** min.
2. holding point at **600°C** holding time: **45** min.
3. holding point at **950°C** holding time: **60** min.

Afterwards the muffel is cooled down to **450°C** in the ARCA 20/30 preheating furnace and is cast after 30 min.

Preheating Biotan Vest C+B

$5^{\circ}\text{C}/\text{min} > 300^{\circ}\text{C}$, 30 min. $5^{\circ}\text{C}/\text{min} > 884^{\circ}\text{C}$, 30 min. $5^{\circ}\text{C}/\text{min} > 450^{\circ}\text{C}$

Re-use

Biotan™ Nb can only be cast once. Do not re-use the alloy.

Crucibles

The copper crucible included with the **DOR-A-MATIC** is used for casting in this system. When casting in the **Biotan casting system**, please use the copper plate.

Casting

The alloy is cast with conventional settings for pure titanium (settings according to the quantities) in the **DOR-A-MATIC** as well as in the **Biotan titanium casting system**.

Afterwards the muffel is cooled at room temperature until the devesting temperature is reached.

Blasting

The cast object is blasted with pure aluminium oxide 125 µm at a maximum of 2 bar and at an obtuse angle with a distance of 5 cm from the object. We recommend tempering the surface with the blasting medium Alpha OKM (Schütz Dental).

Processing

After blasting the surface is worked with the set of special titanium instruments. However, the rotating instruments should only be used for **Biotan™ Nb** then.

Finishing

Finishing and polishing is done with the recommended instruments of the special titanium working set and with sandpaper wheels (Schütz Dental).

Cleaning

After polishing, **Biotan™ Nb must settle for 10 min.** before coming into contact with water or a steamcleaner, in order to allow for a passivation layer to build.

Note!

Biotan™ Nb can also be used for the ceramic facing technique, it was developed e.g. for Tizian™ and VITA Titanium Ceramic.

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SINERGIA DISK Ti are disks made in a Titanium based alloy of excellent quality (Ti6Al4V-ELI) for carrying out dental prostheses by CAD/CAM milling centers.

Ti6Al4V-ELI is a Titanium Grade 5 with highest purity ("Extra-Low Interstitial") and is used since decades in the biomedical field thanks to its excellent biocompatibility, corrosion resistance and physical-mechanical properties.

INDICATIONS

Suitable for abutments, single crowns and long-span bridges on natural dies or implants:

- Anterior: not longer than 3-units between the pillars, connection section > 6 mm²
- Posterior: not longer than 3-units between the pillars, connection section > 9 mm²

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Finishing: use carbide disks (or other proper instruments) to separate the milled framework from the disk. For the finishing use only suitable burs in tungsten carbide with cross-cut and work in the same direction. The milled frameworks can be pre-polished with silicone rubbers and polished with normal polishing pastes and soft brushes (goat hair).

NOTE! During the separation from the disk and the finishing phase, do not overheat the metal frame by avoiding too high speed and pressure. Often cool the piece in water. Never use hydrofluoric acid.

Aesthetic finishing: sandblast with aluminum dioxide 125 - 250 µm at 2-3 bar and further clean with steam. Use composite or ceramic suitable with Titanium. The CTE of the ceramic should be compatible with Titanium ($25-500^{\circ}\text{C} = 9,6 \times 10^{-6} \text{K}^{-1}$), the firing temperature should not be higher than 800°C. Follow the instructions for use of the ceramic manufacturer.

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